

GenoSensor Viral Infection Kit iso-B Home Kit, iso-B and iso-S Kit Components

<u>Kit Components Included:</u>	<u>Qty</u>	<u>Storage</u>
Viral Infection Kit iso-B Home Kit (Catalog# 3301)		
--Sample A, B	400 μ L (24 students) for each	
--2X isothermal amplification Master Mix		Room Temp.
--Staining Solution		Room Temp.
Viral Infection Kit iso-B (Catalog# 3302)		
--Sample A, B	60 μ L for each DNA	-20 $^{\circ}$ C
--2X isothermal amplification Master Mix	270 μ L	-20 $^{\circ}$ C
-- Sample C, D	60 μ L for each DNA	-20 $^{\circ}$ C
--Staining Solution	60 μ L	-20 $^{\circ}$ C
Viral Infection Kit iso-S (Catalog# 3303)		
--Sample A, B	60 μ L for each DNA	-20 $^{\circ}$ C
--Sample C, D	60 μ L for each DNA	-20 $^{\circ}$ C
--2X iso-S Maser Mix(includes enzymes)	270 μ L	-20 $^{\circ}$ C
--DNA Ladderfor gel electrophoresis	30 μ L	-20 $^{\circ}$ C